

Medical Modalities, your reliable connection for AFT beds, introduces the




HydroAire

GENERATION II

Air-Fluidized Therapy (AFT) Bed System



AFT research shows patients heal at a rate 3 times over conventional Group 2 surfaces*



Your
reliable connection
for AFT beds!


Call today for more
information: 1-800-849-2716
www.medicalmodalities.com

- AFT – the Gold Standard in surfaces for patients with non-healing Stage 3 or 4 pressure ulcers on the trunk or pelvis – is a proven, reliable solution for home care.
- Utilization of AFT results in better outcomes and allows agencies to reduce visits and increase healing rates.
- The HydroAire represents the best in AFT: Better fluidization, quieter, creates less heat – the entire surface fluidizes. It is suitable for mobile homes.
Learn more about the HydroAire AFT at www.auroramfg.com
- Reduces electricity costs by 40% over previous unit.
- Medical Modalities is a consistent provider of AFT to Medicare patients.

*Comparison of air-fluidized therapy with other support surfaces used to treat pressure ulcers in nursing home residents. Dr. Ochs, et al. *Ostomy Wound Management*, Feb. 2005.



HydroAire GENERATION II



AFT: The best surface for non-healing Stage 3 or 4 pressure ulcers

- Pay For Performance! AFT is a solution for chronic pressure wounds.
- AFT is best solution for patients who have not improved on Group 2 surfaces, and a potential solution for patients 'capping out' on Group 2 surfaces.
- Patients who have had myocutaneous skin flap procedures are great candidates!
- Reduces risk of acquired ulcers.
- Reduces incidence of re-hospitalization.*
- MM consultants visit with patients before and after delivery to ensure compliance
- MM provides same-day response for servicing beds.



Patient Lift



Trapeze

*The cost of home air-fluidized therapy for pressure sores. A randomized controlled trial. Strauss, et al. *Journal of Family Practice*, July 1991.



Call today for more information: 1-800-849-2716
www.medicalmodalities.com

Product Features for Model 3000

- Adjustable height enhances patient mobility
- Modular design for easy in-home installation
- Breathable GoreTex® sheets
- Dual blower design provides redundancy in event of blower failure
- Adjustable Fluidization
- 84" long x 39" wide; Height adjustable 26" -30"
- Variable heating temp 80° to 100°
- Remote Control to change parameters
- Bed Weight - empty 282 lbs.; w/ beads 1,122 lbs.
- Siderails are standard for safety
- Trapeze bar and stand available
- Works in conjunction with standard patient lifter
- Slant wedge available at varying angles
- Accessible blower assembly allows easier servicing

Specifications

- Maximum recommended patient weight limit 250 lbs.
- Bed Height 30 or 26 inches
- Bed Width
Overall 36 inches
Inside Tank 33-3/4 inches
- Bed Length
Overall 84 inches
Inside Tank 81-3/4 inches
- Bed Weight
Without Micro Spheres 282 lbs.
Micro Spheres 840 lbs.
With Micro Spheres 1,122 lbs.
Max Patient Size 250 lbs.
- Electrical
Voltage 115 VAC
Frequency 60 Hertz
Noise Level <50 DBa
Power Cord Length 10 feet

NOTE: Federal law restricts this device to sale or rental by or on the order of a physician